



**384 Standard Volume - Assay - PCR plates**

 384 BioRad HSP-3801 LV LP Clr Rd Well Rd Btm\*\*\*  
 384 BioRad HSR-4805 LV LP Wh Clr Rd Well Rd Btm\*\*\*

\*\*\*Compounds Australia can only utilise Hard Shell Full Skirted PCR plates - please check with Compounds Australia on labware compatibility

Type	Can Compounds Australia Supply
Assay	No
Assay	No

**384 Low Volume - Storage/Predilution**

 384 Corning 3672 PP LV Sqr Well Rd Btm  
 384 Labcyte LDV Clr Diamond Well  
 384 MatriCal MCR10112 PP White Sqr Well Flt Btm LVol  
 384 Packard 6008890 PP LV Sqr Well Rd Btm

Type	Can Compounds Australia Supply
Storage	Yes
Storage	Yes
Storage	No
Storage	No

**384 Standard Volume - Storage/Predilution**

 384 Axygen 1205QC PP Clr Sqr Well Diamond  
 384 Greiner 781201 PP Clr Sqr Well Flt Btm  
 384 Greiner 784201 PP Small Vol Deep Well U Btm  
 384 Labcyte PPL Clr Sqr Well Flt Btm  
 384 Corning 3656 PP Clr Sqr Well Rd Btm  
 384 Corning 3657 PP Clr Sqr Well Rd Btm  
 384 EP 951031003 PP DW Clr Sqr Well V Btm  
 384 Eppendorf PP 0030624300 V-Btm  
 384 Greiner 781209 PP Blk Sqr Well Flt Btm  
 384 Greiner 781280 PP Sqr Well V Btm  
 384 Matrix 4309 PP Clr Sqr Well V Btm  
 384 Nunc 264573\_264574 PP Clr Sqr Well Flt\_Internal U-Btm  
 384 REMP 180030 PP Clr Sqr Well Diamond

Type	Can Compounds Australia Supply
Storage	Yes
Storage	Yes
Storage	Yes
Storage	Yes
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No
Storage	No

**384 Deep Well - Storage/Predilution**

384 Greiner 781270 PP DW Clr Sqr Well V Btm

Type	Can Compounds Australia Supply
Storage	Yes